

MX-RL-Pro


LCD Digital Tube Rotator, Horizontal

Features

- LCD displays speed and time
- Adjustable speed range of 10 to 70 rpm
- Provides gentle but effective rolling and mixing
- Can be used in cold rooms or incubators
- Wide variety of Disk accessories and Rotisserie accessories available



Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.











Disk Accessories Type-1

Cat. No.18900160

Used for MX-RD-Pro for 1.5mLx60 centrifuge tubes holder



Disk Accessories Type-2 Cat. No.18900161

Used for MX-RD-Pro for 15mLx16 centrifuge tubes holder



Disk Accessories Type-3

Cat. No.18900162

Used for MX-RD-Pro for 50mLx8 centrifuge tubes holder



Disk support rods

Cat. No.18900140

4pcs in a group use with disk accessories







Rotisserie Accessory Type-3

Cat. No.18900144

Used for MX-RL-Pro, for 50mLx24 centrifuge tubes held horizontally



Rotisserie Accessory Type-4

Cat. No.18900145

Used for MX-RL-Pro,for 1.5mLx32 centrifuge tubes held vertically



Rotisserie Accessory Type-1 Cat. No.18900142

Used for MX-RL-Pro, for 1.5mLx48 centrifuge tubes held horizontally



Rotisserie Accessory Type-5 Cat. No.18900146

Used for MX-RL-Pro,for 15mLx16 centrifuge tubes held vertically



Rotisserie Accessory Type-2

Cat. No.18900143

Used for MX-RL-Pro, for 15mLx24 centrifuge tubes held horizontally



Rotisserie Accessory Type-6 Cat. No.18900147

Used for MX-RL-Pro, for 50mLx16 centrifuge tubes held vertically

