



Mini Vortex Mixer

MDB-05F is a hand-held thermostatic dry bath controlled by a microcomputer, and it is featured as its most advantage, super fast heating up. The instruments with small size and light weight are specially designed for use in crowded experimental bench.

GenFollower Biotech CO., LTD.

Address: No.273, Yuexiuzhong road, Cao'e street, Shangyu area, Shaoxing, Zhejiang, China
 Phone: +86-575-82936953 / Fax: +86-575-82032762
 Email: GenFollower@GenFollower.com



⊕ Simple operation

> Push the tube into the soft attachment to start shaking, and stops when tube lifted.

⊕ A variety of containers available

> maximum 50ml centrifuge tube

⊕ Adjustable speed

> maximum speed can reach 4500rpm

⊕ Portable

> Support car cigarette lighter to convert power supply

⊕ Small size, compact and lightweight

> one hand can hold the device

⊕ Two modes are available for choosing

> Pulsing mode for short time mixing sample (inching)
 > Continuous mode (long time continuous operation)

> Technical Parameters

Model	MI-45
Control mode	Pulsing mode or rotation speed adjustable
Amplitude	4.5 mm
Rotating speed	Up to 4500rpm
Maximum support	50ml centrifuge tube
Size (L x W x H)	102mm x 102mm x 88mm
Power supply	DC 12V 1A, 12W; Support car cigarette lighter to convert power supply
Motor Power	1.8W
Weight	0.8kg



We open to new ideas!