



Dry Bath incubator

The GMB-0215 digital dry bath heaters are outstanding products with broad temperature range and accurate control, which also have a patent in preventing from scald. The digital panel of our products is easy and quick to operate, which gives life science researchers the convenience and the efficiency they need when working with precious samples.

GenFollower Biotech CO., LTD.

Address: No.273, Yuexiuzhong road, cao'e street, Shangyu area, Shaoxing, Zhejiang, China
 Phone: +86-575-82936953 / Fax: +86-575-82032762
 Email: GenFollower@GenFollower.com



⊕ **Have a patent in anti-scald**

> When the block temperature reach to 150°C, the temperature of the panel and the body are below 50°C.

⊕ **Rapid heating time**

> Finish heating within 15min

⊕ **Digital touchable keyboard & Full screen**

⊕ **3 steps of programs preset**

> Preset programs can be changed and combined as you like.

⊕ **Wide temperature range**

> Temperature can reach to 150°C

> Technical Parameters

Product Number	GMB-0215
Operation mode	Digital touchable keyboard
Screen	Color LCD screen
Preset program	3 steps of programs selected
Panel temp	<45°C @RT
Temp range	Amb +5° to 150°C (Amb @25°C)
Timer	1s-99h99min99s
Temp accuracy	± 0.3°C @100°C / ± 0.5°C @150°C
Temp uniformity	≤0.4°C
Heating time	≤15min
Heating Power	400W
Temp display accuracy	0.1°C
Dimensions (L x W x H)	202mm x 265.5mm x 86mm
Number of blocks	2
Weight	4.2kg
Voltage/Frequency	220v/50Hz
Warranty	1 Year
Fuse	250V/3A Φ5×20
Plug Type	AC-03



We support your research!