


Shaker Series

Multi-function Shaker (Horizontal+Pendulum+3D)




SYC-FP200 (Horizontal)

Servicebio Multi-functional Shaker



More Prominent Even Mixing Effect

Traditional 3D Shaker



Ineffective Mixing of Middle Components

>> Product Features

- Can be used independently in pendulum mode, horizontal mode, and 3D mode.
- When in 3D mode, you can adjust blending strength individually. Experience true three-dimensional blending for optimal mixing results.
- Automatic leveling function. Timer & power-off memory function.
- Using a stepper motor, the use of stepping motor, precise speed control, and long service life.
- Large tray & digital display screen. The bottom plate is rust-proof treated for greater durability.

Horizontal Shaker (Big Tray)



DS-2H200

>> Product Features

- Timer & Power-off memory functions. Large digital display screen.
- Large tray, large digital display screen; using a stepper motor, precise speed control, and long service life.
- Swinging radius: 15mm.
- The bottom plate is rust-proof treated for greater durability

Pendulum Shaker (Big Tray)



DS-2S100

>> Product Features

- Using a stepper motor, the use of stepping motor, precise speed control, and long service life.
- 10 degrees swing radiation.
- Automatic leveling function. Timing & power-off memory function. Digital display large screen.
- The bottom plate is rust-proof treated for greater durability.

Wuhan Servicebio Technology Co., Ltd.

5th Floor, 22 Building, Biopark, Gaoxin Second Road, 388, East Lake High-Tech Zone, Wuhan, Hubei, China
 027-5111-3188 | www.servicebio.com | info@servicebio.com

Horizontal / Pendulum Shaker (Split-type)



SYC-FP200 (Horizontal)



SYC-FZ100 (Pendulum)

>> Product Features

- Adopting a split-type setting, the flat wire connects the controller, and can be placed in a refrigerator, incubator, chromatography cabinet.
- Internal waterproof treatment, can adapt to high humidity, low temperature and other environments.
- Compact design, using a stepper motor, low noise, more reliable quality.
- Timer function, Power-off memory function.
- Pendulum shaker (SYC-FZ100) has the function of automatic leveling.

Horizontal / Pendulum Shaker



SYC-P200 (Horizontal)



SYC-Z100 (Pendulum)

>> Product Features

- Compact design, using a stepper motor, low noise, more reliable quality.
- Timer function, Power-off memory function.
- Pendulum shaker (SYC-Z100) has the function of automatic leveling.

Horizontal/Pendulum Shaker (Dual-Layer)



SYC-SP200



SYC-SZ100

>> Product Features

- Dual-layer design to accommodate more samples for simultaneous shaking and mixing.
- Multiple adjustable levels of interlayer spacing for placing samples of different heights.
- Automatic leveling to prevent sample spillage; Automatic leveling to prevent sample spillage; power-off memory function, automatically resumes operation after power is restored.
- Equipped with a stepper motor for precise speed control and long service life.

Horizontal/Pendulum Shaker (Clamps)



SYC-WP200



SYC-WZ100

>> Product Features

- Spring grid design, providing more stable sample placement.
- Dual-layer spring bracket allows for holding shorter vessels like culture dishes and taller vessels like conical flasks.
- Automatic leveling to prevent sample spillage; power-off memory function, automatically resumes operation after power is restored.
- Equipped with a stepper motor for precise speed control and long service life.