


Product Specification

| | |
|-------------------|--|
| Item Number | FS-E30 |
| Product Name | PCR Film Sealer (Portable) |
| Dimensions(L*W*H) | 47.2x11.9x14.3cm |
| Size/ctn | 54*18*22cm |
| Weight | 7kg |
| Compatibility | 96-well/384-well PCR reaction plate, ELISA plate. |
| Product Features | <ul style="list-style-type: none"> • Simple operation and reliable sealing • Automatic rolling mode, no need for personnel to operate • One time complete sealing within 30 seconds, high efficiency • Intelligent touch of entrance and exit control; Preventing problems such as samples not being placed properly or personnel making mistakes in operation, timely touch the switch can withdraw the samples from the working chamber • Small footprint, light weight, easy placing and moving at will • The status light can visually judge the machine status. <p>Yellow light : standby state, Green light - running state, Red light - abnormal state, Whenever red light is on, the machine can be restarted directly. And it can be used normally after automatic initialization.</p> |
| Storage Condition | Room Temperature |
| Brand | Servicebio |
| Image |  |

