




### Product Specification

Item Number	FS-A20
Product Name	PCR Film Sealer (Smart)
Dimensions(L*W*H)	36.2*17.8*29.1cm
Size/ctn	43*24*37cm
Weight	10.35kg
Compatibility	96-well/384-well PCR reaction plate, ELISA plate.
Product Features	<ul style="list-style-type: none"> <li>• Touch screen display for easy and intuitive operation.</li> <li>• The sealing film is tight, no evaporation, and good consistency.</li> <li>• Fast and efficient, 20 seconds to complete a flat sealing film.</li> <li>• The parameters are flexible and adjustable, the number of times of rolling (1-9999) and the value of rolling pressure (1-3) can be programmed.</li> <li>• Noise-free operation, creating a comfortable working environment.</li> <li>• The adapter is easy to replace and can be customized according to the experimental requirements.</li> </ul>
Storage Condition	Room Temperature
Brand	Servicebio



3188

