

## Certification of Analysis

**Product Name:** CD3-APC

**Clone:** HIT3a

**Batch Number:** 202209054

**Quantity:** 100T

**Production date:** 20220926

**Expiration Date:** 20240326

**Buffer and Preservative:** 0.01M PBS (pH 7.2) with 0.1% proclin300 and 0.2% BSA

**Method of Test:** Direct immunofluorescence staining for flow cytometric analysis

**Suggested working titer:** 20 $\mu$ l/1 $\times$ 10<sup>6</sup> cells or 100 $\mu$ l peripheral blood

**Results:** This batch of products reacts with normal human peripheral blood lymphocytes (77.4%), is negative in peripheral blood monocytes and granulocytes.

**Conclusion:** This batch of products conforms to the specificity of anti-human CD3 antibody.

Technician: Meng yuan Han Date: 09/26/2022

Authorized for release: Ying Shu Date: 09/26/2022

[www.quantobio.com](http://www.quantobio.com)

E-mail: [info@quantobio.com](mailto:info@quantobio.com)

Tel: +86 10 59773 899

Fax: +86 10 59773 998

*For research use only. Not for diagnostic use. Not for resale. Quantobio will not be held responsible for patent infringement or other violations that may occur with the use of our products.*