

## Wiggins Analytics CON600



### Professional Bench Top Series

The bench-top meters offer high precision measurement for laboratory standards. All meters come with standard electrode stand, probe and electrode and PC link software.

### Professional Bench Top conductivity Meter

- > Large LCD with simultaneous display of conductivity and temperature.
- > Up to 5 calibration points ( at 1 point per range).
- > 99 data memory with Real time clock & date.
- > Computer link with RS232 software kit.

## Specifications

Model No.	CON600 Professional Conductivity Bench meter				
Order Code	BP3020				
Range	0- 19.99 $\mu$ S	0 - 199.9 $\mu$ S	0 -1999 $\mu$ S	0 -19.99 mS	0 - 199.9 mS
Display Resolution	0.01 $\mu$ S	0.1 $\mu$ S	1 $\mu$ S	0.01 mS	0.1 mS
Accuracy	$\pm$ 1% F.S. + 1digit				
Calibration	5 point with 1 point on each range				
TDS factor	Adjustable 0.3 to 1.00				
Temperature compensation	Automatic 0 to 80°C				
Memory	99 data with real time clock				
Operating Temperature	5 to 40°C				
Storage Temperature	-20 to 60°C				
Operating humidity	Up to 95% relative humidity				
Meter Size	217 x 168 x 585mm				
Weight	2.1 kg				



ООО «Компания «ПанЭко»  
115477, город Москва, улица Кантемировская, дом 58, этаж 1 комн 1010  
Телефон: 8-800-550-72-31  
info@paneco-ltd.ru  
www.paneco-ltd.ru