

## Chemical Resistant Diaphragm Pump

### Fast pumping speed for industrial applications


C2000 series chemical resistant diaphragm pumps are ideally suited for pumping aggressive gases and vapors in a vacuum range down to 1mbar and pumping speed up to 245L/min.

#### Features

- > All wetted parts are made of chemical resistant materials
- > Fast pumping speed
- > Low noise
- > Overheat protection
- > Vacuum can be controlled by manual valve, vacuum controller
- > Easy and convenient for maintenance



The motor drives the pump head to vacuum



Model / Specifications	C2000E	C2000T	C2000V
Order No.	W1032001	W1032002	W1032003
Flow rate [m <sup>3</sup> /h] at atm. pressure	14.7	10.8	9.0
Flow rate [L/min] at atm. pressure	245	180	150
Ultimate vacuum [mbar abs.]	< 70	< 2	< 1
Max. Power P [W]	750	750	750
Motor speed [rpm]	1380	1380	1380
Pump head	Single stage	Triple stage	Quadruple stage
Hose connections of inlet	KF25 <sup>1)</sup>	KF25 <sup>1)</sup>	KF25 <sup>1)</sup>
Hose connections of outlet	G1/2 <sup>1)</sup>	G1/2 <sup>1)</sup>	G1/2 <sup>1)</sup>
Noise [dB]	60	60	60
Dimensions W x D x H [mm]	615×285×386	615×285×386	615×285×386
Weight [kg]	41	41	41
Power supply	220~240V/50Hz	220~240V/50Hz	220~240V/50Hz

1) Included : 2 barbed fittings for tubing 16 mm inner dia.