

Remotely Controllable Ultra High Torque / High Speed Stirrers WB3000-EC



- > Suitable for reaction systems and other high torque / high speed applications
- > Set speed can be maintained when viscosity of liquid changes
- > Brushless DC motor made in Germany, for high performance stirring processes
- > Maintenance-free
- > Quiet and reliable
- > Digitally adjustable rotation speed and direction
- > Remote controller can display actual speed and actual torque as well as the set speed
- > Digital (RS-232/485) and analog communication available for remote PC or PLC control

Performance

- > Years of vibration-free and silent operation.
- > Powerful torque capable of handling high viscosity tasks.
- > Maintenance of constant motor speed by PID feedback control system even under conditions of changing viscosity.
- > Best effort function intelligently manages its stirring speed to keep stirring even workload is out of its capacity.
- > Prevention of accidental spills or splashes thanks to microprocessor controlled smooth start and stop functions.



Practical external controller.

- > Intuitive and easy control with touch TFT.
- > It enables convenient and safe external control without opening the sash of the fume hood or safety cabinet. Compact and slim head design for diverse flexibility in configuring other test equipment and accessories. (optional)

Safety

- > Sturdy aluminum main body efficiently absorbing and emitting the heat generated by the motor.
- > Separated adapter from the main body minimizes the risks of electrical hazards to the users.

Specifications

Model	WB3000-EC
Control Type	TFT Display for Speed and Torque and On-Touch Control
Speed Range (rpm)	70-3000
Speed Accuracy(rpm)	±1
Maximum Viscosity (cps)	200000
Maximum Torque (N-cm)	410
Rated Torque (N-cm)	116
Maximum Capacity (L) H2O	100
Power (W)	450
Dimensions (W \times L \times H in mm)	165×75×75
Motor Weight (Kg)	2.8
Order No.	100402

WIGGE	:NS (Gmb	Н	
-------	-------	-----	---	--

Gässlesweg 22-24, 75334 Straubenhardt, Germany

Tel.: 0049 7248 4529088

WIGGENS China

Room 303, Hall C, Office Building M8, No.1 Jiuxiangiao East Road,

Chaoyang District, Beijing 100015, China

Tel: +86 400-809-2068

Fax: +86 400-809-2068-112

Email: info@ wiggens.com service@wiggens.com

Website: www.wiggens.com